

ZOLL

IVTM™ Intravascular Temperature Management



Case Study into Clinical Advantage of IVTM Fever Management using TGXP over Surface Cooling with Covid19 Patient

HOSPITAL SIZE

348 beds

PATIENT HISTORY

Age 67

Weight 94kg / 207lbs

Gender Male

RR 30's

O₂ Sat 80%

Peep 12

FiO₂ 60%

FEVER MANAGEMENT CONSIDERATIONS FOR COVID-19 PATIENTS

ZOLL Temperature Management System Utilized	Thermogard XP® (TGXP)
ZOLL IVTM Catheter Utilized	Quattro® Catheter
Starting Temperature	40.2°C (104.5°F)
Target Temperature	36.5°C (97.7°F)
Time to Target Temperature	<60 minutes (See Figure 1)

IVTM Performance Summary


Patient presented with high fever and respiratory distress syndrome related to COVID-19. Poor fever control is associated with increased mortality in COVID-19 patients.¹

Surface cooling was initiated with no change in patient's temperature and continued deterioration in the patient's overall condition.

While preparing to prone the patient a Quattro® catheter was placed. Patient's hemodynamics stabilized once the **fever was controlled and reduced to 36.5°C – eliminating the need for proning.** Respiratory rate improved to 22, oxygen saturation to 97%, peep to 7, FiO₂ to 30%.

Patient was cooled quickly and reached target temperature in less than an hour. IVTM is a precise and efficient way to control fever.²

Figure 1: IVTM Patient Data



IVTM System Thermogard XP®

Track patient & system data, transferring digitally to patients file post treatment.

Important data are clearly displayed on the large screen

Set upper and lower patient temperature alarm limits

Set the target temperature between 31°C and 38°C

Toggle between Standby and Run modes

Visual alarm indicator



Clearly visible coolant temperature indicator

Control mode, target temperature, and rate setting all on one screen.

Choose a cooling/warming rate



Precise, fast and effective control of patient temperature, reducing clinical time & costs while achieving better results compared to surface cooling.²⁻⁷



Available from



Call 01782 637009 for system details & ordering

¹Tharakan, et al., *Critical Care*. 2020;24:298

²Hoedemaekers CW, et al. *Critical Care*. 2007;11:R91.

³Mayer SA, et al. *Critical Care Medicine*. 2004;32(12):2508-2515.

⁴Diringer MN, et al. *Critical Care Medicine*. 2004;32(12):559-564.

⁵Heard KJ, et al. *Resuscitation*. 2010;81:9-14.

⁶Horn CM, et al. *Journal of Neurointerventional Surgery*. 2014 Mar;6(2):91-95.

⁷Knapik P, et al. *Kardiologia Polska*. 2011;69(11):1157-1163.

