

DORO®

QUICK-CLAMP WITH STARBURST

- FOR RADIOLUCENT QUICK-RAIL® -



BLACK
FOREST
MEDICAL
GROUP

Available from

DELTA

Call 01782 637009 for
system details & ordering

DORO® QUICK-CLAMP WITH STARBURST

- FOR RADIOCULENT QUICK-RAIL® -

The DORO® Quick-Clamp with Starburst for Radiolucent Quick-Rail® serves as an intuitive interface* for devices used in frameless stereotactic or endoscopic procedures.

Designed for use over the sterile barrier and formfittingly connected to the proven Quick-Rail® of the DORO® Skull Clamp Radiolucent, the adaptor provides a stable and secure connection point for third party devices.



Please find further information on the product's web page.



DIRECTLY CONNECT stereotactic devices to the Quick-Rail®

*The interface is designed to bear a max. load of 16.8 Nm (which equals a weight of 3.4 kg attached to a horizontally mounted, 50 cm long arm). It is validated for use with the Brainlab VarioGuide™. It is the responsibility of the user to verify that third party devices can safely be secured to the DORO® Quick-Clamp with Starburst for Radiolucent Quick-Rail®. Product is shown with Brainlab VarioGuide™ and DORO® Skull Clamp Radiolucent.