


DORO® QR3 SKULL CLAMP



Available from
DELTA

Call 01782 637009 for
system details & ordering

DORO® QR3 Skull Clamp



EASY TO READ
laser etched scale
(N and lbs) on the
torque screw.



NAVIGATION READY
Stainless steel starbursts
on both sides of the Quick-
Rail® allow direct attachment
of movement-sensitive
navigation tracking device
adaptors.



QUICK-RAIL®
An easy to locate, easy to handle
interface for a stable and secure
connection of **DORO LUNA®**
and **DORO COBRA®** Retractor
Systems.



**ULTRA-PRECISE
FABRICATION**


No play of two pin rocker arm when locked.

**LARGE DORO® SKULL PIN
PORTFOLIO**

Disposable and reusable skull pins for adult and pediatric applications are available.

**ULTRA-PRECISE
FABRICATION**

Minimum play between skull clamp and extension arm.



DORO® QR3 Skull Clamp
item no. 1001.001



MIX & MATCH

metal and radiolucent parts by means of a **DORO®** Interface Adaptor, item no. 3033-51.

The **DORO® QR3** Skull Clamp provides 3-pin rigid cranial fixation and features the latest generation of Quick-Rail®.

Designed with straight and upright proportions, the QR3 offers more Quick-Rail® space than all its predecessors.

Furthermore, the improved Quick-Rail® design makes it easier to connect accessories, even through multiple layers of sterile drape.



Contact us to find out more and to place an order

☎ 01782 637009

✉ sales@deltasurgical.co.uk

🌐 www.deltasurgical.co.uk

Delta Surgical Ltd, Unit 10, Evolution, Lymedale Business Park, Newcastle-under-Lyme, Staffordshire, ST5 9QF