

Single Use Cranial Access Kit



Cost effective and efficient. Each kit is supplied as a sterile pack

Kit Type

Order Code

UOM

Standard Kit
(detailed overleaf)

RMT0082

10

Standard Kit + CH5 Drill Burr

RMT0082-CH5

10

Cranial Hand Drill

Order Code

UOM

Plastic Cranial Hand Drill

D-5001-S

1

Drills are compatible with cranial access burrs from all manufacturers



DELTA

Delta Surgical Ltd. ☎ 01782 637009 ✉ sales@deltasurgical.co.uk
Unit 10, Evolution, Lymedale Business Park, Newcastle-under-Lyme,
Staffordshire, ST5 9QF



Standard Kit Contents

This convenient, pre-packaged sterile kit provides access to the subarachnoid space. Designed to be used with external drainage and monitoring to reduce intracranial pressure. Each kit contains the following:

- Disposable Hand Drill
- Bone Wax
- Suture Mersilk 3/8 Circle RC
- 26mm, black 2/0, 45cm
- Suture
- Skin Marker Pen, Fine (with Ruler)
- Syringe, 10ml-12ml, LS, Eccentric
- Bard Parker Handle with No 11 Blade
- Bard Parker Handle with No 15 Blade
- Holzheimer Blunt
- Mosquito Forceps, Curved
- Mayo Needle Holder, 16cm
- Adson Tissue Forceps 12cm
- Adson Dressing Forceps
- Iris Scissor, 11cm
- Blade Scalpel Carbon Steel No 10
- Gruenwald Bayonet Forceps, Serrated Tip, 8inch
- MacDonald Dissector, Curved, 185 mm (7 1/4"), Double Ended, Blunt/Blunt
- Fenestrated Drape 60 x 50cm
- Sponge Stick
- Absorbent Towel
- Drape Trolley/Table Cover 112x112

