




# THE NEXT FRONTIER IN THROMBECTOMY

The first and only  
integrated Micro-Filtration  
technology




designed to  
maximize Clot Retention

# NeVa™

**DROP ZONE™ TECHNOLOGY**  
DESIGNED FOR 1ST PASS SUCCESS  
WITH ALL CLOT TYPES

# NeVa™ NET

**THE FIRST AND ONLY**  
MICRO FILTRATION TECHNOLOGY  
DESIGNED TO MAXIMIZE CLOT RETENTION



Product No	Description	Diameter (mm)	Working Length (mm)	Total Basket Length (mm)	Pusher-wire Length (cm)	Number of Drop Zones	Minimum MC ID (inches)
30020V-MS	NeVa 4.0 x 22 mm	40	22	39	180	2	0.021
VN-4529-03RR	NeVa 4.5 x 29 mm	45	29	46	180	3	0.021
VN-5537-03RR	NeVa 5.5 x 37 mm	55	37	56	180	3	0.027
VN-4030-03NR	NeVa NET 4.0 x 30 mm	40	30	49	180	3	0.027
VN-5537-03NR	NeVa NET 5.5 x 37 mm	55	37	56	180	3	0.027