



NEUROMEDEX<sup>®</sup>

VentriFix<sup>®</sup>  
Fixation for Ventricular Catheters

## VentriFix®

### Innovative Fixation for Ventricular Catheters

VentriFix® is a novel product, designed to fix ventricular drains securely.

VentriFix® largely eliminates the risk of catheter dislocation by unintended tensile stress, occurring, for example, during post-operative transfer or lifting of the patient.

VentriFix® specific design reduces micromovements of the ventricular drainage catheter – an essential contribution to infection prophylaxis and wound healing.

VentriFix® enhances the prevention of catheter-associated infections and thus application of our ventricular catheters VentriGuard® and VentriClear®.

## VentriFix® fixation for ventricular catheters

- Unique product design for catheter fixation
- Reduced infection risk during external ventricular drainage
- Increased efficiency: reduced risk of revision surgery due to ventricular drain dislocation
- Ideal complement to the VentriGuard® and VentriClear® product concepts



### VentriFix® Fixation for Ventricular Catheters

Ref. No.	Description	Size	Packaging Unit
59-0030	VentriFix® Fixation for Ventricular Catheters	8.5 F	10 ea/box
available from the 4th quarter 2019	VentriFix® Fixation for Ventricular Catheters	10 F / 12 F	10 ea/box

### VentriGuard® Ventricular Catheter

Ref. No.	Description	Size	I.D. / O.D.	Length	Perforation	Packaging Unit
60-0200	VentriGuard® Ventricular Catheter	8.5 F	1.6 mm / 2.8 mm	30 cm	standard	1 ea/box
60-0010	VentriGuard® Ventricular Catheter	10 F	1.9 mm / 3.3 mm	30 cm	large	1 ea/box

### VentriClear® Ventricular Catheter

Ref. No.	Description	Size	I.D. / O.D.	Length	Perforation	Packaging Unit
60-1210	VentriClear® Ventricular Catheter	8.5 F	1.6 mm / 2.8 mm	30 cm	standard	10 ea/box
60-1010	VentriClear® Ventricular Catheter	10 F	1.9 mm / 3.3 mm	30 cm	large	10 ea/box
60-1510	VentriClear® Ventricular Catheter	12 F	2.3 mm / 4.0 mm	30 cm	large	10 ea/box

Length markings in 1 cm intervals



Visual puncture control due to semi-transparent catheter material







Neuromedex GmbH    phone +49 (0)40 696 564 101  
Vierenkamp 15      fax    +49 (0)40 696 564 200  
D 22453 Hamburg    web    www.neuromedex.com