

## L-Line MEMORY Microinstruments





1. | MICRO-WORKING ENDS

- A variety of styles available for different surgical approaches
- Delicate working ends for precise performance


4. | 360° ROTATABLE WORKING ENDS

- Hand controlled wheel ensures exact positioning of the working end
- Rotation combined with malleable shaft offers unparalleled instrument positioning





6. | COLOUR CODING

- Simple identification of the working ends by colour coated handles



2. | MEMORY SHAFT (NITINOL)

- Shapes to coincide with anatomical needs
- Reverts to original shape after sterilization



3. | RINSING-CONNECTION/CLEANING

- Built in flushing port for easy and efficient cleaning of internal components

5. | ERGONOMIC HANDLE

- Balanced and low weight handle for fatigue-free working



# L-Line MEMORY Microinstruments



**46.600** (90 mm)

**46.610** (110 mm)

Micro-Scissors, straight,  
sharp, movable one side,  
Colour: yellow



**46.620** (90 mm)

**46.630** (110 mm)

Micro-Scissors, sharp,  
straight, Colour: yellow



**46.601** (90 mm)

**46.611** (110 mm)

Micro-Scissors, sharp,  
curved, Colour: yellow



**46.623** (90 mm)

**46.633** (110 mm)

Micro-Scissors, blunt,  
straight, movable one side,  
Colour: yellow



**46.622** (90 mm)

**46.632** (110 mm)

Micro-Scissors, blunt,  
straight, Colour: yellow



**46.621** (90 mm)

**46.631** (110 mm)

Micro-Scissors, blunt,  
curved, Colour: yellow



**46.625** (90 mm)

**46.635** (110 mm)

Micro-Scissors, sharp,  
45° angled, movable  
one side, Colour: yellow



**46.604** (90 mm)

**46.614** (110 mm)

Micro-Needleholder, fine,  
curved, Colour: blue



**46.605** (90 mm)

**46.615** (110 mm)

Micro-Needleholder, raw,  
curved, Colour: blue



**46.609** (90 mm)

**46.619** (110 mm)

Micro-Forcep, serrated,  
Ø 3 mm, Colour: black



**46.602** (90 mm)

**46.612** (110 mm)

Micro-Forcep, tip 0.5 mm,  
Colour: green



**46.603** (90 mm)

**46.613** (110 mm)

Micro-Forcep, tip 0.9 mm,  
Colour: green



**46.606** (90 mm)

**46.616** (110 mm)

Micro-Spoon, sharp,  
Ø 2 mm, Colour: red



**46.607** (90 mm)

**46.617** (110 mm)

Micro-Spoon, sharp,  
Ø 3 mm, Colour: red



**46.608** (90 mm)

**46.618** (110 mm)

Micro-Spoon, sharp,  
Ø 5 mm, Colour: red



**46.626** (90 mm)

**46.636** (110 mm)

Blakesley, sharp,  
Ø 3 mm, Colour: red