

Manufacture and sales:

VOSTRA GmbH      Phone +49 241 96850-0  
medical products      Fax +49 241 96850-56  
Im Süsterfeld 3      info@vostra.de  
52072 Aachen      www.vostra.de  
Germany

CE 0481




Vostra GmbH was founded in 1985 with the goal of developing, manufacturing and marketing high-quality consumables for the neurosurgery, ENT and general surgery.

Our customers needs have always been the focus in our strive for excellence. Dedication and expertise have nourished an international customer base which we supply on short notice – world wide.

Vostra – excellence in medical products.

## Spherical Neuro Swabs

Round swabs for neurosurgery



# Spherical Neuro Swabs


## Product Features

### Round swabs for neurosurgery

Spherical neuro swabs have been developed for applications in neurosurgery and are manufactured in a specialized production process that guarantees the product to retain its shape during use.

Round neuro swabs are supplied in standard form with a reinforced thread of surgical silk. The sewing procedure used during the manufacturing process places the end of the thread in the interior of the product, guaranteeing maximum product safety. The tear-resistant attachment of the thread and the x-ray contrast material inside the spherical neuro swab make it simple and safe to use. Spherical neuro swabs are also optionally available without thread.

Spherical neuro swabs have a wide range of uses and should be a standard part of any neurosurgical kit. The surgeon can tailor the product to his/her needs and determine the x-ray contrast marking, size and set quantities.



### Features

#### Safety

The attached thread and the sewn-in x-ray marking ensure a safe application.

#### Material properties

Complex manufacturing processes guarantee the material to retain its shape.

#### Customised range of sizes

Products and designs are manufactured according to customers specifications.

#### Quality

Spherical neuro swabs are manufactured under meticulously defined cleanroom conditions. Maximum product safety is guaranteed via continuous microbiological process control.

# Spherical Neuro Swabs

## Sizes



Spherical neuro swabs with x-ray markings, double-sterile pouch

Spherical neuro swabs with thread, double-sterile

Sets of 5

Size	Article no.
10 mm	10950510
15 mm	10950515
20 mm	10950520
25 mm	10950525

Sets of 10

Size	Article no.
10 mm	10951010
15 mm	10951015
20 mm	10951020
25 mm	10951025

Packaging unit

30 sets per box

Spherical neuro swabs, double-sterile

Sets of 5

Size	Article no.
10 mm	10970510
15 mm	10970515
20 mm	10970520
25 mm	10970525

Sets of 10

Size	Article no.
10 mm	10971010
15 mm	10971015
20 mm	10971020
25 mm	10971025

