

SKATING ON THIN ICE?

ZOLL IVTM™ Intravascular Temperature Management

Thermogard XP® Key Facts

- ✓ Helps you avoid “afterdrop,” an abrupt drop in core temperature from the redistribution of heat to peripheral tissues
- ✓ Maintains target temperature within $\pm 0.2^{\circ}\text{C}$ 97% of the time
- ✓ Offers cost-effective, multifunction catheters with built-in CVC functionality
- ✓ Track patient & system data, transferring digitally to patients file post treatment



Available from



Call 01782 637009 for system details & ordering

ZOLL IVTM System Thermogard XP®

Precise, fast and effective control of patient temperature, reducing clinical time & costs

Important data are clearly displayed on the large screen

Set upper and lower patient temperature alarm limits

Set the target temperature between 31°C and 38°C



Toggle between Standby and Run modes

Visual alarm indicator

Clearly visible coolant temperature indicator

Control mode, target temperature, and rate setting all on one screen.

Choose a cooling/warming rate



Available from
DELTA

Call 01782 637009 for system details & ordering