Single-Use

Cranial Access Kit

CE


Kit Type	Order Code	UOM
Standard Kit	RMT0082-DS	10
Standard Kit + CH5 Drill Burr* *Standard Kit + CH5 Drill Burr is not available for export	RMT0082-CH5-DS	10

Additional Sterile, Single-Use Skull ICP Drill Burrs with adjustable collar and key	Order Code	UOM
2.7mm ICP metal skull bolt drill burr	S-AI-DBS-2700	1
4.5mm Raumedic CH5 ICP metal skull bolt drill burr	S-AI-DBS-4500	1
5.8mm plastic ICP skull bolt & EVD drill burr	S-AI-DBS-5800	1
6.5mm Raumedic CH9 ICP metal skull bolt drill burr	S-AI-DBS-6500	1

- Bone Wax LN903
- Suture LN1625 24 nylus blu 45cm m2 usp3 single armed ps
- Suture LN321 silk blk 45cm m3 usp2 0 single arm fs
- Swab Gauze 10x10 XRD 12ply tied 5
- Gallipot 120ml
- Skin Marker Pen Fine (with ruler)
- Syringe, 10ml, LS, Eccentric Emerald
- Scalpel with No 11 Blade
- Scalpel with No 15 Blade
- Razor for Hair
- BD Blunt Fill Needle with Filter Red 18G x 1 1/2"
- MacDonald Dissector Curved 185mm double ended blunt/blunt
- Sponge Stick 20cm Orange
- Drape Opaque Film 100x100cm with S/A & Adhesive Fenstration 10cm
- Disposable Hand Drill with variable chuck
- Holzheimer Blunt
- Mosquito Forceps, Curved
- Mayo Needle Holder, 16cm
- Adson Tissue Forceps 12cm
- Adson Dressing Forceps
- Iris Scissors, 11cm
- Blade CS No 10
- Gruenwald Bayonet Forceps, serrated tip, 8inch length
- Drape 75 x 75cm Fen 6x8cm S/A SMS
- BD Microlance Orange Needle 25G x 5/8"
- Spinal Needle 18g x 3.5"
- Absorbent Towel

Cost effective and efficient. Each kit is supplied as a sterile pack



A selection of the kit items listed are displayed here

